## ZnS(Ag) Zinc Sulfide Scintillation Material

Silver activated zinc sulfide has a very high scintillation efficiency comparable to that of NaI(TI). It is only available as a polycrystalline powder.

ZnS(Ag) has a maximum in the scintillation emission spectrum at 450nm. The light conversion efficiency is relatively poor for fast electrons which may be an advantage when detecting heavy ionizing particles in a relatively intense  $\gamma$ -ray background. Scintillation decay times between several hundreds of ns and 10µs are reported and phosphorescence of still longer duration has been noted.

Thicknesses greater than about 25 mg/cm<sup>2</sup> become unusable because of the opacity of the multicrystalline layer to its own luminescence.

Its use is limited to thin screens used primarily for  $\alpha$ -particles or other heavy ion detection. The main application of ZnS(Ag) is in  $\alpha$ -particle monitors.

The figure below shows the energy spectrum response for 5.5 MeV

incident  $\alpha$ -particles. The broad peak is clearly above the electronic noise.

ZnS(Ag) can be used to detect thermal neutrons if a lithium compound enriched in <sup>6</sup>Li is incorporated. The alpha particle and triton from the <sup>6</sup>Li (n,  $\alpha$ ) <sup>3</sup>H reaction produce scintillations upon interacting with the ZnS(Ag). The figure also shows a thermal neutron spectrum.

Another use for ZnS(Ag) is for the detection of fast neutrons. A fast neutron detector is made by imbedding the ZnS(Ag) in a clear hydrogenous compound. The neutrons are detected by measuring the recoiling proton from a neutron proton scattering interaction. The protons will have a broad spectrum of energies up to the energy of the incident neutrons.



## Figure 1. Spectra for thermal neutron and 5.5 MeV alpha response of ZnS(Ag)

# **BICRON**<sup>®</sup>

## Properties –

Density [g/cm³] 4.09
Cleavage plane polycrystalline
Wavelength of emission max. [nm] 450
Lower wavelength cutoff [nm]
Refractive index @ emission max 2.36
Photoelectron yield [% of Nal(Tl)] (for γ-rays) 130
Decay constant [ns] 110

SAINT-GOBAIN CRYSTALS & DETECTORS **Scintillation Products** 

## ZnS(Ag) Zinc Sulfide Scintillation Material

SAINT-GOBAIN CRYSTALS & DETECTORS

## USA

Saint-Gobain Crystals & Detectors 12345 Kinsman Road Newbury, OH 44065 Tel: (440) 564-2251 Fax: (440) 564-8047

#### Europe

**Saint-Gobain Cristaux & Detecteurs** 104 Route de Larchant BP 521 77794 Nemours Cedex, France Tel: 33 (1) 64 45 10 10 Fax: 33 (1) 64 45 10 01

## Japan

**Saint-Gobain Crystals & Detectors KK** 3-7, Kojimachi, Chiyoda-ku, Tokyo 102-0083 Japan Tel: 81 (0) 3 3263 0559 Fax: 81 (0) 3 5212 2196

## China

**Saint-Gobain China Investment Co., Ltd.** 24-05 CITIC Building 19 Jianguomenwai Ave. Beijing 100004 China Tel: 86 (0) 10 6513 0311 Fax: 86 (0) 10 6512 9843

www.detectors.saint-gobain.com